

For Research Use Only

CoraLite® Plus 488-conjugated G6PD Recombinant monoclonal antibody

Catalog Number: CL488-86240



Basic Information

Catalog Number: CL488-86240	GenBank Accession Number: BC000337	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2539	CloneNo.: 250831F8
Isotype: IgG	UNIPROT ID: P11413	Recommended Dilutions: IF/ICC: 1:50-1:500
Immunogen Catalog Number: AG21862	Full Name: glucose-6-phosphate dehydrogenase	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
	Calculated MW: 59 kDa	
	Observed MW: 60 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : NIH/3T3 cells,
Species Specificity: human, mouse	

Background Information

G6PD (Glucose-6-phosphate 1-dehydrogenase) belongs to the glucose-6-phosphate dehydrogenase family. It is a housekeeping enzyme encoded in mammals by an X-linked gene. G6PD has important functions in intermediary metabolism because it catalyzes the first step in the pentose phosphate pathway and provides reductive potential in the form of NADPH. (PMID:7489710). This protein has 3 isoforms produced by alternative splicing.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed NIH/3T3 cells using CoraLite® Plus 488 G6PD antibody (CL488-86240, Clone: 250831F8) at dilution of 1:200.