

For Research Use Only

FcZero-rAb® CoraLite® Plus 488 Anti-Mouse CD80 Rabbit Recombinant Antibody

Catalog Number: **CL488-FcA98158**



Basic Information

Catalog Number: CL488-FcA98158	GenBank Accession Number: NM_009855.2	Purification Method: Protein A purification
Concentration: 100ug, 500 ug/ml	GeneID (NCBI): 12519	CloneNo.: 241371E12
Source: Rabbit	UNIPROT ID: Q00609-1	Recommended Dilutions: FC: 0.25 ug per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: CD80 antigen	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: EG0840	Calculated MW: 35 kDa	

Applications

Tested Applications: FC	Positive Controls: FC : mouse peritoneal macrophages,
Species Specificity: mouse	

Background Information

CD80 (also known as B7-1) is a type I membrane protein that is a member of the immunoglobulin superfamily, with an extracellular immunoglobulin constant-like domain and a variable-like domain required for receptor binding. It is expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, monocytes, and macrophages. CD80 is the receptor for the proteins CD28 and CTLA-4 found on the surface of T-cells. It is involved in the costimulatory signal essential for T-lymphocyte activation. T-cell proliferation and cytokine production is induced by the binding of CD28, binding to CTLA-4 has opposite effects and inhibits T-cell activation. CD80 also acts as a cellular attachment receptor for adenovirus subgroup B. (PMID: 7545666; 12015893; 16920215)

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse peritoneal macrophages were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse CD80 Rabbit RecAb (CL488-FcA98158, Clone:241371E12), or 0.25 ug FcZero-rAb® CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10⁶ mouse peritoneal macrophages were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse CD80 Rabbit RecAb (CL488-FcA98158, Clone:241371E12)(red), or 0.25 ug FcZero-rAb® CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-FcA98136, Clone: 240953C9)(blue). Cells were not fixed.