

For Research Use Only

FcZero-rAb® CoraLite® Plus 488 Anti-Mouse PD-L1/CD274 Rabbit Recombinant Antibody



Catalog Number: CL488-FcA98166-4

Basic Information

Catalog Number:
CL488-FcA98166-4

Concentration:
100ug, 500 ug/ml

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
NM_021893.3

GeneID (NCBI):
60533

UNIPROT ID:
Q9EP73

Full Name:
CD274 antigen

Calculated MW:
33 kDa

Purification Method:
Protein A purification

CloneNo.:
242265E6

Recommended Dilutions:
FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:
493 nm / 522 nm

Applications

Tested Applications:
FC

Species Specificity:
mouse

Positive Controls:

FC : anti-CD3/CD28 treated mouse splenocytes,

Background Information

Programmed cell death ligand 1 (PD-L1, CD274, or B7-H1) is the first member of B7 family to be discovered. B7 family molecules are type I transmembrane proteins belonging to the immunoglobulin superfamily. In concert with their CD28 family receptors, the B7s are key regulators of the adaptive immune response. PD-L1 is suggested as a negative regulator of T and B cells, and plays an important role in mediating tolerance of lymphocytes to self-antigens. It is also involved in the costimulatory signal, essential for T-cell proliferation and production of IL-10 and IFN gamma, in an IL-2-dependent and a PD-1-independent manner.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or anti-CD3/CD28 treated mouse splenocytes were surface stained with CoraLite® Plus 647 Anti-Mouse CD3e (145-2C11) and 0.25 ug CoraLite® Plus 488 Anti-Mouse PD-L1/B7-H1 Rabbit RecAb (CL488-FcA98166-4, Clone: 242265E6). Cells were not fixed.

1x10⁶ anti-CD3/CD28 treated mouse splenocytes were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse PD-L1/B7-H1 Rabbit RecAb (CL488-FcA98166-4, Clone: 242265E6) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (CL488-FcA98136, Clone: 240953C9) (blue). 1x10⁶ untreated mouse splenocytes were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse PD-L1/B7-H1 Rabbit RecAb (CL488-FcA98166-4, Clone: 242265E6) (black, dashed). Cells were not