

For Research Use Only

# CoraLite®555-conjugated TFAP2C Polyclonal antibody



Catalog Number:CL555-14572

## Basic Information

Catalog Number:

CL555-14572

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5648

GenBank Accession Number:

BC035664

GeneID (NCBI):

7022

UNIPROT ID:

Q92754

Full Name:

transcription factor AP-2 gamma  
(activating enhancer binding protein  
2 gamma)

Calculated MW:

49 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima  
wavelengths:

557 nm / 570 nm

## Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> MCF-7 cells were intracellularly stained with 0.4 ug CoraLite® 555 Anti-Human TFAP2C (CL555-14572) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).