

For Research Use Only

CoraLite®555-conjugated TNF Alpha Monoclonal antibody



Catalog Number:CL555-69002

Basic Information

Catalog Number:

CL555-69002

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1014

GenBank Accession Number:

GeneID (NCBI):

7124

ENSEMBL Gene ID:

ENSG00000232810

Full Name:

tumor necrosis factor (TNF
superfamily, member 2)

Purification Method:

Protein G purification

CloneNo.:

11H8A5

Excitation/Emission maxima
wavelengths:

557 nm / 570nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human

Background Information

TNF, as also known as TNF-alpha, or cachectin, is a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. It is expressed as a 26 kDa membrane bound protein and is then cleaved by TNF-alpha converting enzyme (TACE) to release the soluble 17 kDa monomer, which forms homotrimers in circulation. It is produced chiefly by activated macrophages, although it can be produced by many other cell types such as CD4+ lymphocytes, NK cells, neutrophils, mast cells, eosinophils, and neurons. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, INS resistance, and cancer.

This antibody can be used to neutralize the bioactivity of TNF alpha.

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ PMA, LPS and Brefeldin A treated THP-1 cells were intracellularly stained with 0.25 ug CoraLite®555 Anti-Human NeutraKine® TNF Alpha (CL555-69002, Clone:11H8A5) (red), or 0.25 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Intracellular Staining Permeabilization Wash Buffer.