

For Research Use Only

CoraLite® Plus 555-conjugated NeuN Recombinant antibody

Catalog Number: CL555-84401-4



Basic Information

Catalog Number:

CL555-84401-4

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25689

GenBank Accession Number:

NM_001082575

GeneID (NCBI):

146713

UNIPROT ID:

A6NFN3

Full Name:

hexaribonucleotide binding protein 3

Observed MW:

46-52 kDa

Purification Method:

Protein A purification

CloneNo.:

241461G7

Recommended Dilutions:

IF-P: 1:50-1:500

Excitation/Emission maxima

wavelengths:

554 nm / 570 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human, mouse, rat, pig

Positive Controls:

IF-P : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse brain tissue using CoraLite® Plus 555 NeuN antibody (CL555-84401-4, Clone: 241461G7) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).