

For Research Use Only

CoraLite®594-conjugated OSBP Polyclonal antibody



Catalog Number:CL594-11096

Featured Product

Basic Information

Catalog Number:

CL594-11096

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1578

GenBank Accession Number:

BC011581

GeneID (NCBI):

5007

UNIPROT ID:

P22059

Full Name:

oxysterol binding protein

Calculated MW:

89 kDa

Observed MW:

89 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug Coralite®594 Anti-Human OSBP (CL594-11096) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).