

For Research Use Only

CoraLite®594-conjugated RBP7 Polyclonal antibody



Catalog Number: **CL594-14541**

Basic Information

Catalog Number:

CL594-14541

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6045

GenBank Accession Number:

BC063013

GeneID (NCBI):

116362

UNIPROT ID:

Q96R05

Full Name:

retinol binding protein 7, cellular

Calculated MW:

15.5 kDa

Observed MW:

15-18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human, mouse, rat

Positive Controls:

IF-P : human kidney tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human kidney tissue using CoraLite® 594 RBP7 antibody (CL594-14541) at dilution of 1:200.