

For Research Use Only

CoraLite®594-conjugated FSHR Polyclonal antibody



Catalog Number: CL594-22665 **3 Publications**

Basic Information

Catalog Number:

CL594-22665

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18158

GenBank Accession Number:

BC118548

GeneID (NCBI):

2492

UNIPROT ID:

P23945

Full Name:

follicle stimulating hormone receptor

Calculated MW:

695 aa, 78 kDa

Observed MW:

71-78 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Cited Applications:

IF, IHC

Species Specificity:

human, mouse

Cited Species:

duck

Positive Controls:

IF : HeLa cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jia Qi	33129993	Transl Res	IHC
Yu Zhang	34222406	Front Vet Sci	IF
Yu Zhang	34421865	Front Microbiol	IF

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CL594-22665 (FSHR antibody) at dilution of 1:100.