

For Research Use Only

CoraLite®594-conjugated YTHDF2 Polyclonal antibody



Catalog Number:CL594-24744

Featured Product

Basic Information

Catalog Number:

CL594-24744

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5922

GenBank Accession Number:

BC002559

GeneID (NCBI):

51441

UNIPROT ID:

Q9Y5A9

Full Name:

YTH domain family, member 2

Calculated MW:

62 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

Species Specificity:

human, mouse, rat

Background Information

YTHDF2 also named as CLL associated antigen KW 14 or HGRG8 is a 579 amino acid protein, which contains one YTH domain and exists as two alternatively spliced isoforms. YTHDF2 is expressed in in pancreas, testis and placenta. YTHDF2 may play a role in human longevity.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

