

For Research Use Only

CoraLite®594-conjugated Cadherin-17 Monoclonal antibody



Catalog Number:CL594-60351

Basic Information

Catalog Number:

CL594-60351

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19488

GenBank Accession Number:

BC113464

GeneID (NCBI):

1015

UNIPROT ID:

Q12864

Full Name:

cadherin 17, LI cadherin (liver-intestine)

Calculated MW:

832 aa, 92 kDa

Observed MW:

120 kDa

Purification Method:

Protein G purification

CloneNo.:

3B10A9

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human, pig

Positive Controls:

IF-P : human colon cancer tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human colon cancer tissue using CoraLite®594 Cadherin-17 antibody (CL594-60351, Clone: 3B10A9) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human colon cancer tissue using CoraLite®594 Cadherin-17 antibody (CL594-60351, Clone: 3B10A9) at dilution of 1:200.