

For Research Use Only

CoraLite®594 Anti-Mouse LPAM-1 (DATK32)

Catalog Number:CL594-65101



Basic Information

Catalog Number:

CL594-65101

Concentration:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC011184

GeneID (NCBI):

16421

Full Name:

integrin beta 7

Purification Method:

Affinity purification

CloneNo.:

DATK32

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse bone marrow cells were surface stained with 0.5 ug CoraLite® 594 Anti-Mouse LPAM-1 (CL594-65101, Clone:DATK32) or 0.5 ug Isotype Control. Cells were not fixed.