

For Research Use Only

# CoraLite®594 Anti-Mouse LPAM-1 (DATK32)



Catalog Number:CL594-65101

## Basic Information

Catalog Number:

CL594-65101

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC011184

GeneID (NCBI):

16421

Full Name:

integrin beta 7

Purification Method:

Affinity purification

CloneNo.:

DATK32

Excitation/Emission maxima  
wavelengths:

588 nm / 604 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse bone marrow cells were surface stained with 0.5 ug CoraLite® 594 Anti-Mouse LPAM-1 (CL594-65101, Clone:DATK32) or 0.5 ug Isotype Control. Cells were not fixed.