

For Research Use Only

CoraLite®594 Anti-Mouse CD41 (MWReg30)



Catalog Number: **CL594-65189**

Basic Information

Catalog Number:

CL594-65189

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

BC120493

GeneID (NCBI):

16399

UNIPROT ID:

Q9QUM0

Full Name:

integrin alpha 2b

Purification Method:

Affinity purification

CloneNo.:

MWReg30

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 Balb/c mouse peripheral blood platelets were surface stained with 0.5 ug CoraLite@594 Anti-Mouse CD41 (CL594-65189, Clone: MWRReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.



1×10^6 Balb/c mouse peripheral blood platelets were surface stained with 0.5 ug CoraLite@594 Anti-Mouse CD41 (CL594-65189, Clone: MWRReg30), or 0.5 ug Isotype Control. Cells were not fixed.