

For Research Use Only

CoraLite®594 Anti-Rat CD86 (24F)

Catalog Number: **CL594-65224**



Basic Information

Catalog Number:

CL594-65224

Size:

100ug, 500 µg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

NM_020081

GeneID (NCBI):

56822

UNIPROT ID:

O35531

Full Name:

CD86 molecule

Calculated MW:

36 kDa

Purification Method:

Affinity purification

CloneNo.:

24F

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

Rat

Background Information

CD86 (also known as B7.2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and cytotoxic T-lymphocyte-associated protein 4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with cytotoxic T-lymphocyte-associated protein 4 negatively regulates T-cell activation and diminishes the immune response.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ unstimulated or LPS-stimulated (5 ug/mL, 3 days) rat splenocytes were surface co-stained with 0.5 ug CoraLite® Plus 647 Anti-Rat CD45RA (OX-33) and 0.5 ug CoraLite® 594 Anti-Rat CD86 (CL594-65224, Clone: 24F). Cells were not fixed.