

For Research Use Only

# CoraLite® Plus 594 Mouse IgG2a Isotype Control (MOPC-173)



Catalog Number: **CL594-65244**

## Basic Information

Catalog Number:

CL594-65244

Size:

100ug, 0.2 mg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

N/A

CloneNo.:

MOPC-173

Excitation/Emission maxima  
wavelengths:

594 nm / 615 nm

## Applications

Tested Applications:

FC

Species Specificity:

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite® Plus 594 Anti-Mouse IgG2a Isotype Control (CL594-65244, Clone:MOPC-173) (red), unstained (blue), or CL594 Anti-Human CD3 (Hit3a) (green). Cells were not fixed. Lymphocytes were gated.