

For Research Use Only

CoraLite®594 Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: CL594-65571



Basic Information

Catalog Number:

CL594-65571

Concentration:

100ug, 500 ug/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

n/a

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were unstained (blue) or surface stained with 5 ul CoraLite®594 Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody (CL594-65571, Clone: MOPC-173) (red) or CoraLite® Plus 594 Anti-Human CD3 (Hit3a) (CL594-65112, Clone: Hit3a) (green). Cells were not fixed.