

CoraLite®594-conjugated ERp57/ERp60 Monoclonal antibody

Catalog Number: **CL594-66423**

Basic Information

Catalog Number:

CL594-66423

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8741

GenBank Accession Number:

BC014433

GeneID (NCBI):

2923

UNIPROT ID:

P30101

Full Name:

protein disulfide isomerase family A, member 3

Calculated MW:

505 aa, 57 kDa

Observed MW:

57 kDa, 60 kDa

Purification Method:

Protein G purification

CloneNo.:

3G4G7

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima

wavelengths:
588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : U2OS cells,

Background Information

PDIA3, also named as P58, ER60, ERp57, ERp60, ERp61, GRP57, GRP58 and PI-PLC, is a member of the PDI family, participates in the oxidation, reduction, and isomerization of disulfide bonds for correct folding of secretory proteins before modification and transport in the endoplasmic reticulum. It is associated with apoptosis or inhibition of cancer cell growth. PDIA3 was once thought to be a phospholipase; however, it has been demonstrated that this protein actually has protein disulfide isomerase activity. It is thought that complexes of lectins and PDIA3 mediate protein folding by promoting formation of disulfide bonds in their glycoprotein substrates.

Storage

Storage:


Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using CoraLite® 594 ERp57/ERp60 antibody (CL594-66423, Clone: 3G4G7) at dilution of 1:200.