

For Research Use Only

CoraLite®594-conjugated HRPT2, CDC73 Monoclonal antibody



Catalog Number:CL594-66490

Featured Product

Basic Information

Catalog Number:

CL594-66490

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2959

GenBank Accession Number:

BC014351

GeneID (NCBI):

79577

UNIPROT ID:

Q6P1J9

Full Name:

cell division cycle 73, Paf1/RNA
polymerase II complex component,
homolog (S. cerevisiae)

Calculated MW:

531 aa, 61 kDa

Observed MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

4B4H2

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human, Pig, Mouse, Rat

Positive Controls:

IF-P : human breast cancer tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human breast cancer tissue using CL594-66490 (HRPT2; CDC73 antibody) at dilution of 1:50.



Immunofluorescent analysis of (4% PFA) fixed human breast cancer tissue using CL594-66490 (HRPT2; CDC73 antibody) at dilution of 1:50.