

For Research Use Only

CoraLite®594-conjugated CHCHD5 Monoclonal antibody



Catalog Number:CL594-66584

Basic Information

Catalog Number:

CL594-66584

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22534

GenBank Accession Number:

BC004498

GeneID (NCBI):

84269

UNIPROT ID:

Q9BSY4

Full Name:

coiled-coil-helix-coiled-coil-helix
domain containing 5

Calculated MW:

110 aa, 12 kDa

Observed MW:

14 kDa

Purification Method:

Protein G purification

CloneNo.:

2G10A7

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human

Positive Controls:

IF/ICC : MCF-7 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite®@594 CHCHD5 antibody (CL594-66584, Clone: 2G10A7) at dilution of 1:200.