

For Research Use Only

CoraLite®594-conjugated RAB3GAP1 Monoclonal antibody

Catalog Number: CL594-66688

Featured Product



Basic Information

Catalog Number:

CL594-66688

Size:

1000 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17710

GenBank Accession Number:

BC146809

GeneID (NCBI):

22930

UNIPROT ID:

Q15042

Full Name:

RAB3 GTPase activating protein subunit 1 (catalytic)

Calculated MW:

981 aa, 111 kDa

Observed MW:

130 kDa

Purification Method:

Protein G purification

CloneNo.:

4A3A1

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC, FC (Intra)

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : HEK-293 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using Coralite®594 RAB3GAP1 antibody (CL594-66688, Clone: 4A3A1) at dilution of 1:200.



1X10⁶ HEK-293 cells were intracellularly stained with 0.4 ug Coralite®594 Anti-Human RAB3GAP1 (CL594-66688, Clone:4A3A1) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).