

For Research Use Only

CoraLite®594-conjugated USP4 Monoclonal antibody



Catalog Number:CL594-66822

Basic Information

Catalog Number:

CL594-66822

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28297

GenBank Accession Number:

BC125131

GeneID (NCBI):

7375

UNIPROT ID:

Q13107

Full Name:

ubiquitin specific peptidase 4 (proto-oncogene)

Calculated MW:

963 aa, 109 kDa

Observed MW:

109 kDa

Purification Method:

Protein A purification

CloneNo.:

1G11G8

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

Species Specificity:

Human, Mouse

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data