

For Research Use Only

CoraLite®594-conjugated PRTN3 Monoclonal antibody



Catalog Number:CL594-67030

Basic Information

Catalog Number:

CL594-67030

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG25915

GenBank Accession Number:

BC096183

GeneID (NCBI):

5657

UNIPROT ID:

P24158

Full Name:

proteinase 3

Calculated MW:

256 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

2D1B8

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima

wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

Human

Positive Controls:

IF-P : human tonsillitis tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 PRTN3 antibody (CL594-67030, Clone: 2D1B8) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 PRTN3 antibody (CL594-67030, Clone: 2D1B8) at dilution of 1:200.