

For Research Use Only

# CoraLite®594-conjugated RALB Monoclonal antibody



Catalog Number:CL594-67094

## Basic Information

Catalog Number:

CL594-67094

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28600

GenBank Accession Number:

BC018163

GeneID (NCBI):

5899

UNIPROT ID:

P11234

Full Name:

v-ral simian leukemia viral oncogene homolog B (ras related; GTP binding protein)

Calculated MW:

206 aa, 23 kDa

Observed MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

1D11C11

Recommended Dilutions:

IF 1:200-1:800

Excitation/Emission maxima

wavelengths:  
588 nm / 604 nm

## Applications

Tested Applications:

IF-P

Species Specificity:

Human, rat, pig, mouse

Positive Controls:

IF : human lung cancer tissue,

## Background Information

## Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human lung cancer tissue using CoraLite@594 RALB antibody (CL594-67094, Clone: 1D11C11 ) at dilution of 1:400.