

For Research Use Only

CoraLite®594-conjugated RAB3D Monoclonal antibody



Catalog Number:CL594-67118

Basic Information

Catalog Number: CL594-67118	GenBank Accession Number: BC016471	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 9545	CloneNo.: 3C2F10
Source: Mouse	UNIPROT ID: O95716	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2a	Full Name: RAB3D, member RAS oncogene family	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG16784	Calculated MW: 219 aa, 24 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF : A549 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using CoraLite®594 RAB3D antibody (CL594-67118, Clone: 3C2F10) at dilution of 1:200.