

For Research Use Only

CoraLite®594-conjugated ARL2BP Monoclonal antibody



Catalog Number:CL594-67316

Basic Information

Catalog Number:

CL594-67316

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29261

GenBank Accession Number:

BC003087

GeneID (NCBI):

23568

UNIPROT ID:

Q9Y2Y0

Full Name:

ADP-ribosylation factor-like 2 binding protein

Calculated MW:

19 kDa

Observed MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1G5C2

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite®594 Anti-Human ARL2BP (CL594-67316, Clone:1G5C2) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).