

For Research Use Only

CoraLite®594-conjugated TTPA Monoclonal antibody



Catalog Number:CL594-67431

Basic Information

Catalog Number: CL594-67431	GenBank Accession Number: BC058000	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 7274	CloneNo.: 3E7D10
Source: Mouse	UNIPROT ID: P49638	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2b	Full Name: tocopherol (alpha) transfer protein	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG13587	Calculated MW: 32 kDa	
	Observed MW: 32 kDa	

Applications

Tested Applications: IF-P	Positive Controls: IF : human liver cancer tissue,
Species Specificity: Human, Mouse, Rat, Pig	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite®594 TTPA antibody (CL594-67431, Clone: 3E7D10) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite®594 TTPA antibody (CL594-67431, Clone: 3E7D10) at dilution of 1:200.