

For Research Use Only

CoraLite®594-conjugated SNF8 Monoclonal antibody



Catalog Number:CL594-67696

Basic Information

Catalog Number:

CL594-67696

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30384

GenBank Accession Number:

BC008976

GeneID (NCBI):

11267

UNIPROT ID:

Q96H20

Full Name:

SNF8, ESCRT-II complex subunit,
homolog (S. cerevisiae)

Calculated MW:

258 aa, 29 kDa

Observed MW:

29 kDa

Purification Method:

Protein G purification

CloneNo.:

2C10F2

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human

Positive Controls:

IF : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CoraLite®594-conjugated SNF8 antibody (CL594-67696, Clone: 2C10F2) at dilution of 1:100.