

For Research Use Only

# CoraLite® Plus 594-conjugated NSFL1C Monoclonal antibody



Catalog Number:CL594-68182

## Basic Information

Catalog Number:

CL594-68182

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8040

GenBank Accession Number:

BC002801

GeneID (NCBI):

55968

UNIPROT ID:

Q9UNZ2

Full Name:

NSFL1 (p97) cofactor (p47)

Calculated MW:

41 kDa

Observed MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1F12E12

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima  
wavelengths:

594 nm / 615 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

Human, Mouse, Rat, Pig

Positive Controls:

IF : A431 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A431 cells using CoraLite® Plus 594 NSFL1C antibody (CL594-68182, Clone: 1F12E12 ) at dilution of 1:200.