

CoraLite®594-conjugated Mono/Di-Methyl-Histone H3 (Lys9) Recombinant antibody

Catalog Number: **CL594-80219**

Basic Information

Catalog Number:

CL594-80219

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066245

GeneID (NCBI):

8350

UNIPROT ID:

P68431

Full Name:

histone cluster 1, H3a

Observed MW:

15-18 kDa

Purification Method:

Protein A purification

CloneNo.:

1E12

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

Species Specificity:

Human, mouse, rat

Background Information

Histones, including H1/H5 (linker histones), H2, H3, and H4 (core histones), are nucleic proteins which interact with DNA to form the nucleosomes and play important roles in gene regulation and DNA replication. Histone proteins are highly post-translationally modified while Histone H3 is the most extensively modified. Methylation of Histone H3 at lysine 9 is linked to transcriptional repression. This antibody is specific to monomethyl or dimethyl-Histone H3 while it does not recognize trimethyl-Histone H3 (Lys9). It is also named as H3K9me1/2.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

