

For Research Use Only



CoraLite®594-conjugated Cardiac Troponin T Recombinant monoclonal antibody

Catalog Number:CL594-81638-3

Basic Information

Catalog Number: CL594-81638-3	GenBank Accession Number: BC002653	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7139	CloneNo.: 25023081
Isotype: IgG	UNIPROT ID: P45379	Recommended Dilutions: IF-P: 1:50-1:500
Immunogen Catalog Number: AG7716	Full Name: troponin T type 2 (cardiac)	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
	Calculated MW: 36 kDa	
	Observed MW: 34-40 kDa	

Applications

Tested Applications: IF-P	Positive Controls: IF-P : mouse heart tissue,
Species Specificity: human, mouse, rat, pig	

Background Information

Cardiac Troponin T, also named as cTnT, TNNT2. This gene encodes the cardiac isoform of troponin T. The encoded protein is the tropomyosin-binding subunit of the troponin complex, which is located on the thin filament of striated muscles and regulates muscle contraction in response to alterations in intracellular calcium ion concentration. Mutations in this gene have been associated with familial hypertrophic cardiomyopathy as well as with dilated cardiomyopathy. Cardiac Troponin T can be detected as about 34-40 kDa.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse heart tissue using CoraLite®594 Cardiac Troponin T antibody (CL594-81638-3, Clone: 250230B1) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse heart tissue using CoraLite®594 Cardiac Troponin T antibody (CL594-81638-3, Clone: 250230B1) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).