

For Research Use Only

CoraLite®594-conjugated Phospho-EGFR (Tyr1197) Recombinant antibody

Catalog Number:CL594-84906



Basic Information

Catalog Number:

CL594-84906

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC094761

GeneID (NCBI):

1956

UNIPROT ID:

P00533

Full Name:

epidermal growth factor receptor
(erythroblastic leukemia viral (v-erb-
b) oncogene homolog, avian)

Calculated MW:

1210 aa, 134 kDa

Observed MW:

175 kDa

Purification Method:

Protein A purification

CloneNo.:

241826G10

Recommended Dilutions:

FC (Intra): 0.25 ug per 10⁶ cells in a
100 µl suspension

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Positive Controls:

FC (Intra) : EGF treated A431 cells,

Background Information

Epidermal Growth Factor Receptor (EGFR) is one of the receptor tyrosine kinases (RTKs), members of the ErbB/HER family which consists of ErbB1 (EGFR or HER1), ErbB2 (HER2 or Neu), ErbB3 (HER3), and ErbB4 (HER4). The EGFR autophosphorylation sites Tyr-992 and Tyr-1173 (also known as Tyr 1197) have been shown to play a critical role in the activation of the MAPK cascade following EGF stimulation. (PMID: 15708852)

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ A431 cells untreated (dashed lines) or treated with EGF which intracellularly stained with 0.25 ug Coralite® 594 Phospho-EGFR (Tyr1197) Recombinant Antibody (CL594-84906, Clone:241826G10)(red), or 0.25 ug Coralite® 594 Rabbit IgG Isotype Control RecAb (CL594-98136, Clone: 240953C9) (blue). Cells were fixed with 4% PFA and permeabilized with 90% MeOH.