

For Research Use Only

# CoraLite®594-conjugated PCNP Recombinant monoclonal antibody

Catalog Number:CL594-86075



## Basic Information

Catalog Number:

CL594-86075

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30951

GenBank Accession Number:

BC013916

GeneID (NCBI):

57092

UNIPROT ID:

Q8WW12

Full Name:

PEST proteolytic signal containing  
nuclear protein

Calculated MW:

19 kDa

Observed MW:

24-28 kDa

Purification Method:

Protein A purification

CloneNo.:

250576B6

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima  
wavelengths:

588 nm / 604 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite®594 PCNP antibody (CL594-86075, Clone: 250576B6) at dilution of 1:200, CL488-phalloidin (green).