

For Research Use Only



CoraLite®594 Anti-Human SIRP Alpha/CD172a Rabbit Recombinant Antibody

Catalog Number:CL594-98101

Basic Information

Catalog Number: CL594-98101	GenBank Accession Number: NM_080792	Purification Method: Protein A purification
Concentration: 100tests, 5 ul/test	GeneID (NCBI): 140885	CloneNo.: 240789E1
Source: Rabbit	UNIPROT ID: P78324	Recommended Dilutions: FC: 5 ul per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: signal-regulatory protein alpha	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: EG0625	Calculated MW: 55 kDa	

Applications

Tested Applications: FC	Positive Controls: FC : human PBMCs,
Species Specificity: human	

Background Information

SIRP Alpha, also known as CD172a, SHPS-1, and BIT, belongs to the SIRP family. SIRP Alpha is a transmembrane glycoprotein expressed explicitly on myeloid cells and provides a "do not eat me" signal after engaging with integrin-associated protein CD47 on tumor cells (PMID: 36419386). SIRP Alpha shows heterogeneity in molecular weight (~65-120 kDa) in various tissues due to differential glycosylation (PMID: 18051954).

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 5 ul CoraLite®594 Anti-Human SIRP Alpha/CD172a Rabbit RecAb (CL594-98101, Clone: 240789E1) or CoraLite®594 Rabbit IgG Isotype Control Recombinant Antibody (CL594-98136, Clone: 240953C9). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



1x10⁶ human PBMCs were surface stained with 5 ul CoraLite®594 Anti-Human SIRP Alpha/CD172a Rabbit RecAb (CL594-98101, Clone: 240789E1) (red) or CoraLite®594 Rabbit IgG Isotype Control Recombinant Antibody (CL594-98136, Clone: 240953C9) (blue). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.