

For Research Use Only

CoraLite® Plus 647-conjugated TFAP2C Polyclonal antibody



Catalog Number: CL647-14572

Featured Product

Basic Information

Catalog Number:

CL647-14572

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5648

GenBank Accession Number:

BC035664

GeneID (NCBI):

7022

UNIPROT ID:

Q92754

Full Name:

transcription factor AP-2 gamma
(activating enhancer binding protein
2 gamma)

Calculated MW:

49 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.2 ug CoraLite® Plus 647 Anti-Human TFAP2C (CL647-14572) (red), or 0.2 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).