

For Research Use Only

# CoraLite® Plus 647-conjugated TRAF1 Polyclonal antibody



Catalog Number:CL647-26845

## Basic Information

Catalog Number:

CL647-26845

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25397

GenBank Accession Number:

BC024145

GeneID (NCBI):

7185

UNIPROT ID:

Q13077

Full Name:

TNF receptor-associated factor 1

Calculated MW:

416 aa, 46 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> Raji cells were intracellularly stained with 0.4 ug CoraLite® Plus 647 Anti-Human TRAF1 (CL647-26845) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).