

For Research Use Only

# CoraLite® Plus 647 Anti-Mouse CD115 (c-fms) (AFS98)



Catalog Number: CL647-65095

## Basic Information

Catalog Number:

CL647-65095

Size:

100 tests, 5 µl/test

Source:

Rat

Isotype:

Rat IgG2a, kappa

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

Excitation/Emission maxima  
wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> Balb/c mouse peritoneal macrophages were surface stained with 5 ul CoraLite® Plus 647-conjugated Anti-Mouse CD115 (c-fms) (CL647-65095, Clone: AFS98) (blue) or CoraLite® Plus 647-conjugated Rat IgG2a isotype control (black). Cells were not fixed.