

For Research Use Only

CoraLite® Plus 647 Anti-Mouse CD24 (M1/69)



Catalog Number: **CL647-65121**

Basic Information

Catalog Number:

CL647-65121

Size:

100ug, 500 µg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with FITC-Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and 0.5 ug CoraLite® Plus 647-conjugated Anti-Mouse CD24 (CL647-65121, Clone: M1/69). Cells were not fixed.



1X10⁶ C57BL/6 mouse splenocytes were surface stained with FITC-Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and CoraLite® Plus 647-conjugated rat IgG2b isotype control antibody. Cells were not fixed.