

For Research Use Only

CoraLite® Plus 647 Anti-Human STRO-1 (STRO-1)



Catalog Number: **CL647-65184**

Basic Information

Catalog Number:

CL647-65184

Size:

100 tests, 5 µl/test

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Purified by precipitation and chromatography

CloneNo.:

STRO-1

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human whole blood cells were surface stained with 5 ul CoraLite® Plus 647 Anti-Human STRO-1 (CL647-65184, Clone: STRO-1). Cells were not fixed. CD45dim, CD34+ cells were gated.