

For Research Use Only

CoraLite® Plus 647 Anti-Mouse CD94 (18D3)



Catalog Number: **CL647-65240**

Basic Information

Catalog Number:

CL647-65240

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor, subfamily D, member 1

Purification Method:

Affinity purification

CloneNo.:

18D3

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 1 ug CoraLite® Plus 647 Anti-Mouse CD94 (CL647-65240, Clone:18D3) (red), or 1 ug Control Antibody. Cells were not fixed.