

For Research Use Only

CoraLite® Plus 647 Anti-Mouse CD49b (HMA2)

Catalog Number: CL647-65245



Basic Information

Catalog Number:

CL647-65245

Size:

100ug, 500 µg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC065139

GeneID (NCBI):

16398

UNIPROT ID:

Q62469

Full Name:

integrin alpha 2

Purification Method:

Affinity purification

CloneNo.:

HMA2

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse NK1.1 (CD161) (FITC-65138, Clone: PK136) and 0.5 ug CoraLite® Plus 647 Anti-Mouse CD49b (CL647-65245, Clone: HMA 2) or 0.5 ug CoraLite® Plus 647-conjugated Armenian Hamster IgG. Cells were not fixed.