

For Research Use Only

# CoraLite® Plus 647 Anti-Mouse CD49b (HMA2)



Catalog Number: **CL647-65245**

## Basic Information

Catalog Number:

CL647-65245

Size:

100ug, 500 µg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC065139

GeneID (NCBI):

16398

UNIPROT ID:

Q62469

Full Name:

integrin alpha 2

Purification Method:

Affinity purification

CloneNo.:

HMA2

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with FITC Anti-Mouse NK1.1 (CD161) (FITC-65138, Clone: PK136) and 0.5 ug CoraLite® Plus 647 Anti-Mouse CD49b (CL647-65245, Clone: HM<sup>Δ</sup> 2) or 0.5 ug CoraLite® Plus 647-conjugated Armenian Hamster IgG. Cells were not fixed.