

For Research Use Only

# CoraLite® Plus 647 Anti-Human CD96 (NK92.39)



Catalog Number: CL647-65250

## Basic Information

Catalog Number:

CL647-65250

Size:

100tests, 5 µl/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC020749

GeneID (NCBI):

10225

UNIPROT ID:

P40200

Full Name:

CD96 molecule

Calculated MW:

402 aa, 46 kDa

Purification Method:

Affinity purification

CloneNo.:

NK92.39

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface co-stained with Anti-Human CD19 and 5 ul CoraLite® Plus 647 Anti-Human CD96 (CL647-65250, Clone: NK92.39), or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.