

For Research Use Only



CoraLite® Plus 647-conjugated CD107b / LAMP2 Recombinant monoclonal antibody

Catalog Number: **CL647-86269-3**

Basic Information

Catalog Number: CL647-86269-3	GenBank Accession Number: NM_002294.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3920	CloneNo.: 250948D7
Isotype: IgG	ENSEMBL Gene ID: ENSG00000005893	Recommended Dilutions: IF/ICC: 1:50-1:500
Immunogen Catalog Number: EG4998	UNIPROT ID: P13473-1	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
	Full Name: lysosomal-associated membrane protein 2	
	Calculated MW: 45 kDa	
	Observed MW: 100-130 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : U-251 cells,
Species Specificity: human	

Background Information

LAMP2 (CD107b) is a Lysosomal membrane glycoprotein. LAMP2 is extensively glycosylated with asparagine-linked oligosaccharides which protect it from intracellular proteolysis (PMID: 10521503). Although LAMP-2 is localized primarily in the endosome-lysosomal membrane of cells, it is also found on the plasma membrane under certain circumstances, e.g., after platelet activation, during granulocytic differentiation and activation, and in some tumor cells (PMID: 12221139). LAMP is involved in lysosomal stability and autophagy (PMID: 12221139). This glycoprotein provides selectins with carbohydrate ligands. LAMP2 may plays a role in tumor cell metastasis (PMID 9426697).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed U-251 cells using CoraLite® Plus 647 CD107b / LAMP2 antibody (CL647-86269-3, Clone: 250948D7) at dilution of 1:200, CL488-phalloidin (green).