

For Research Use Only

FcZero-rAb® CoraLite® Plus 647 Anti-Human IL-4 Rabbit Recombinant Antibody

Catalog Number: CL647-FcA98138



Basic Information

Catalog Number:

CL647-FcA98138

Concentration:

100 tests, 5 µl/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0104

GenBank Accession Number:

BC070123

GeneID (NCBI):

3565

UNIPROT ID:

P05112

Full Name:

interleukin 4

Calculated MW:

153 aa, 17 kDa

Purification Method:

Protein A purification

CloneNo.:

240794E11

Recommended Dilutions:

FC: 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : PMA and ionomycin treated human PBMCs,

Background Information

IL-4 is a cytokine produced by CD4⁺ T cells in response to helminthes and other extracellular parasites. It promotes the proliferation and differentiation of antigen presenting cells. IL-4 also plays a pivotal role in antibody isotype switching and stimulates the production of IgE. This cytokine has been applied in the treatment of autoimmune disorder like multiple myeloma, cancer, psoriasis, and arthritis. IL-4 has also been extensively applied to inhibit detrimental effect of Th1. It may promote the growth of epithelial tumors by mediating increased proliferation and survival.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ PMA and ionomycin treated or untreated human PBMCs were surface stained with FcZero-rAb® CoraLite® Plus 405 Anti-Human CD3 and 5 ul CoraLite® Plus 647 Anti-Human IL-4 Rabbit RecAb (CL647-FcA98138, Clone:240794E11). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).