

FcZero-rAb[®] CoraLite[®] Plus 647 Anti-Human IFN-gamma Rabbit Recombinant Antibody

Catalog Number: **CL647-FcA98187**

Basic Information

Catalog Number:

CL647-FcA98187

Concentration:

100tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0004

GenBank Accession Number:

BC070256

GeneID (NCBI):

3458

ENSEMBL Gene ID:

ENSG00000111537

UNIPROT ID:

P01579

Full Name:

IFN gamma

Calculated MW:

166 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

241390B10

Recommended Dilutions:FC (Intra): 5 ul per 10⁶ cells in a 100 μ l suspension**Excitation/Emission maxima wavelengths:**

654 nm / 674 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Positive Controls:

FC (Intra) : PMA and ionomycin treated human PBMCs,

Background Information

Interferon gamma (IFN γ), is a type II interferon that provides immunity against bacterial, viral and protozoan infections. In its active form, IFN γ is a glycosylated, non-covalently linked homodimer of 29-32 kDa subunits. It is produced by a number of immune cell types including natural killer cells, natural killer T cells, and effector lymphocyte T cells following antigenic and inflammatory triggers. The IFN γ dimer binds to its cognate receptor which has two subunits: IFN- γ R1 which is the ligand-binding chain (α chain) and IFN- γ R2, the signal-transducing chain (β chain). Binding to the receptor activates the JAK/STAT pathway which in turn activates IFN γ responsive genes. While IFN γ can inhibit viral replication, it also works as an immune-modulator and immune-stimulator by increasing surface expression of class I MHC proteins. (PMID: 19268625; 10688427)

Storage


Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

Selected Validation Data



1x10⁶ untreated or PMA and ionomycin treated human PBMCs were intracellularly stained with FcZero-rAb® CoraLite® Plus 488 Anti-Human CD3 (OKT3) and 5 ul CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit RecAb (CL647-FcA98187, Clone: 241390B10). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).

1x10⁶ PMA and ionomycin treated human PBMCs were intracellularly stained with 5 ul CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit RecAb (CL647-FcA98187, Clone: 241390B10) (red) or FcZero-rAb® CoraLite® Plus 647 Rabbit IgG Isotype Control Recombinant Antibody (CL647-FcA98136, Clone: 240953C9) (blue). 1x10⁶ untreated human PBMCs were intracellularly stained with 5 ul CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit RecAb (CL647-FcA98187,

