

For Research Use Only



# CoraLite®647-conjugated Transgelin-2-specific Polyclonal antibody

Catalog Number: ~~CL650-15508~~

Featured Product

## Basic Information

Catalog Number:  
CL650-15508

Size:  
1000 µg/ml

Source:  
Rabbit

Isotype:  
IgG

GenBank Accession Number:  
NM\_003564

GeneID (NCBI):  
8407

UNIPROT ID:  
P37802

Full Name:  
transgelin 2

Calculated MW:  
22 kDa

Observed MW:  
22 kDa

Purification Method:  
Antigen affinity purification

Excitation/Emission maxima  
wavelengths:  
654 nm / 674 nm

## Applications

Tested Applications:  
FC (Intra)

Species Specificity:  
Human

## Background Information

The transgelin family is a group of proteins that belong to 22 kDa actin-related corpnin superfamily. Of all three isoforms, transgelin 1 is the best characterized. Transgelin 1, also known as SM22 alpha, is a specific marker for differentiated smooth muscle cells. Transgelin 2, also known as SM22 beta, is expressed by both smooth muscle and non-smooth muscle cells in a temporally and spatially regulated pattern. Transgelin 3, also known as NP25, is only found in highly differentiated neuronal cells. This antibody can bind transgelin 2 specifically, but not transgelin 1 or 3.

## Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data