

For Research Use Only

CoraLite®647-conjugated Cytokeratin 7 Polyclonal antibody

Catalog Number:CL650-22208



Basic Information

Catalog Number: CL650-22208	GenBank Accession Number: BC002700	Purification Method: Antigen affinity purification
Size: 1000 µg/ml	GeneID (NCBI): 3855	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Source: Rabbit	UNIPROT ID: P08729	
Isotype: IgG	Full Name: keratin 7	
Immunogen Catalog Number: AG17528	Calculated MW: 469 aa, 51 kDa	
	Observed MW: 51 kDa	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
human

Background Information

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are classified into two major sequence types. Type I keratins are a group of acidic intermediate filament proteins, including K9-K23, and the hair keratins Ha1-Ha8. Type II keratins are the basic or neutral counterparts to the acidic type I keratins, including K1-K8, and the hair keratins, Hb1-Hb6. KRT7 is a type II keratin. It is specifically expressed in the simple epithelia lining the cavities of the internal organs and in the gland ducts and blood vessels.

Storage

Storage:
Store at -20°C. Avoid exposure to light.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data