

For Research Use Only

CoraLite®647 Rat IgG2a Isotype Control (2A3)



Catalog Number:CL650-65209

Basic Information

Catalog Number:

CL650-65209

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Affinity purification

Size:

100ug, 200 µg/ml

Full Name:

CloneNo.:

2A3

Source:

Rat

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Isotype:

IgG2a, κ

Applications

Tested Applications:

FC

Species Specificity:

Background Information

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

