

For Research Use Only

CoraLite®647 Mouse IgM Isotype Control (MOPC-104E)



Catalog Number:CL650-65243

Basic Information

Catalog Number:

CL650-65243

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

N/A

Size:

100ug, 200 µg/ml

Full Name:

CloneNo.:

MOPC-104E

Source:

Mouse

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Isotype:

IgM, λ

Applications

Tested Applications:

FC

Species Specificity:

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

